

## ELSPBP1

---

**Reactivity:**Human

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:**26kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human ELSPBP1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

E12; HE12; EDDM12;

**Background:**

The protein encoded by this gene belongs to the sperm-coating protein family of epididymal origin. This protein and its canine homolog are the first known examples of proteins with four tandemly arranged fibronectin type 2 (Fn2) domains in the Fn2-module protein family. [provided by RefSeq, Jul 2008]

*To place an order, please [Click HERE](#).*

**Catalog #:**A2151

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**64100

**Isotype:**IgG

**Swiss Prot:**Q96BH3

**Purity:**Affinity purification

For research use only.