

SCN4B

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:25kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human SCN4B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

Gm1471

Catalog #:A2154

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:6330

Isotype:IgG

Swiss Prot:Q7M729

Purity:Affinity purification

For research use only.

To place an order, please [Click HERE](#).