

## IRF3

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC IF IP ChIP

**Calculated MW:**427aa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic Peptide of human IRF3

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

ks

**Catalog #:**A2172

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**3661

**Isotype:**IgG

**Swiss Prot:**Q14653

**Purity:**Affinity purification

For research use only.

*To place an order, please [Click HERE](#).*