

CUX1

Reactivity:Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:1000 - 1:2000

Calculated MW:77kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human CUX1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

k

Synonym:

CDP; CUX; p75; CASP; CDP1; COY1; Clox; p100; p110; p200; CUTL1; GOLIM6; CDP/Cut; Cux/CDP; Nbla10317;

Catalog #:A2213

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1523

Isotype:IgG

Swiss Prot:P39880

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the homeodomain family of DNA binding proteins. It may regulate gene expression, morphogenesis, and differentiation and it may also play a role in the cell cycle progression. Several alternatively spliced transcript variants encoding different isoforms have been identified.

To place an order, please [Click HERE](#).