

CHD4

Reactivity:Human

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:221kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human CHD4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

Mi-2b; Mi2-BETA;

Catalog #:A2222

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1108

Isotype:IgG

Swiss Prot:Q14839

Purity:Affinity purification

For research use only.

Background:

The product of this gene belongs to the SNF2/RAD54 helicase family. It represents the main component of the nucleosome remodeling and deacetylase complex and plays an important role in epigenetic transcriptional repression. Patients with dermatomyositis develop antibodies against this protein.

To place an order, please [Click HERE](#).