

KAT2A

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:200 - 1:2000

Calculated MW:94kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human KAT2A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

1µg/ml

Synonym:

GCN5; hGCN5; GCN5L2; PCAF-b;

Catalog #:A2224

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:2648

Isotype:IgG

Swiss Prot:Q15758

Purity:Affinity purification

For research use only.

Background:

KAT2A, or GCN5, is a histone acetyltransferase (HAT) that functions primarily as a transcriptional activator. It also functions as a repressor of NF-kappa-B (see MIM 164011) by promoting ubiquitination of the NF-kappa-B subunit RELA (MIM 164014) in a HAT-independent manner (Mao et al., 2009)

To place an order, please [Click HERE](#).