

SUV39H1

Reactivity:Human Mouse Rat

Tested applications:WB IP

Recommended Dilution:WB 1:500 - 1:2000 IP 1:20 - 1:50

Calculated MW:48kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human SUV39H1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A2237

Antibody Type:

Monoclonal Antibody

Species:23

Gene ID:6839

Isotype:IgG

Swiss Prot:O43463

Purity:Affinity purification

For research use only.

Background:

This gene encodes an evolutionarily-conserved protein containing an N-terminal chromodomain and a C-terminal SET domain. The encoded protein is a histone methyltransferase that trimethylates lysine 9 of histone H3, which results in transcriptional gene silencing. Loss of function of this gene disrupts heterochromatin formation and may cause chromosome instability. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2013]

To place an order, please [Click HERE](#).