

## TRIM28

**Reactivity:**Human Mouse Rat

**Tested applications:**WB

**Recommended Dilution:**WB 1:1000 - 1:3000

**Calculated MW:**89kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human TRIM28

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

fl

**Synonym:**

KAP1; TF1B; RNF96; TIF1B;

**Catalog #:**A2245

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**10155

**Isotype:**IgG

**Swiss Prot:**Q13263

**Purity:**Affinity purification

For research use only.

**Background:**

The protein encoded by this gene mediates transcriptional control by interaction with the Kruppel-associated box repression domain found in many transcription factors. The protein localizes to the nucleus and is thought to associate with specific chromatin regions. The protein is a member of the tripartite motif family. This tripartite motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region.

*To place an order, please [Click HERE](#).*