

EZH1

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:200 - 1:1000

Calculated MW:85kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human EZH1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A2258

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:2145

Isotype:IgG

Swiss Prot:Q92800

Purity:Affinity purification

For research use only.

Background:

EZH1 is a component of a noncanonical Polycomb repressive complex-2 (PRC2) that mediates methylation of histone H3 (see MIM 602812) lys27 (H3K27) and functions in the maintenance of embryonic stem cell pluripotency and plasticity (Shen et al., 2008 [PubMed 19026780]).[supplied by OMIM, Mar 2009]

To place an order, please [Click HERE](#).