

MEN1

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100 IF 1:50 - 1:100

Calculated MW: 68kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human MEN1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

MEAI; SCG2;

Catalog #: A2264

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4221

Isotype: IgG

Swiss Prot: O00255

Purity: Affinity purification

For research use only.

Background:

This gene encodes menin, a putative tumor suppressor associated with a syndrome known as multiple endocrine neoplasia type 1. In vitro studies have shown menin is localized to the nucleus, possesses two functional nuclear localization signals, and inhibits transcriptional activation by JunD, however, the function of this protein is not known. Two messages have been detected on northern blots but the larger message has not been characterized. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).