

OGG1

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:200 - 1:1000

Calculated MW:47kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human OGG1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

ef

Synonym:

HMMH; MUTM; OGH1; HOGG1;

Catalog #:A2268

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:4968

Isotype:IgG

Swiss Prot:O15527

Purity:Affinity purification

For research use only.

Background:

This gene encodes the enzyme responsible for the excision of 8-oxoguanine, a mutagenic base byproduct which occurs as a result of exposure to reactive oxygen. The action of this enzyme includes lyase activity for chain cleavage. Alternative splicing of the C-terminal region of this gene classifies splice variants into two major groups, type 1 and type 2, depending on the last exon of the sequence. Type 1 alternative splice variants end with exon 7 and type 2 end with exon 8. All variants share the N-terminal region in common, which contains a mitochondrial targeting signal that is essential for mitochondrial localization. Many alternative splice variants for this gene have been described, but the full-length nature for every variant has not been determined.

To place an order, please [Click HERE](#).