

## SMARCD1

**Reactivity:**Human

**Tested applications:**WB

**Recommended Dilution:**WB 1:200 - 1:2000

**Calculated MW:**58kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic peptide of human SMARCD1

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

f

**Synonym:**

Rsc6p; BAF60A; CRACD1;

**Catalog #:**A2273

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**6602

**Isotype:**IgG

**Swiss Prot:**Q96GM5

**Purity:**Affinity purification

For research use only.

**Background:**

The protein encoded by this gene is a member of the SWI/SNF family of proteins, whose members display helicase and ATPase activities and which are thought to regulate transcription of certain genes by altering the chromatin structure around those genes. The encoded protein is part of the large ATP-dependent chromatin remodeling complex SNF/SWI and has sequence similarity to the yeast Swp73 protein. Two transcript variants encoding different isoforms have been found for this gene.

**To place an order, please [Click HERE](#).**