

## EDF1

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**16kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human EDF1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

fj

**Synonym:**

MBF1; EDF-1;

**Catalog #:**A2283

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**8721

**Isotype:**IgG

**Swiss Prot:**O60869

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes a protein that may regulate endothelial cell differentiation, lipid metabolism, and hormone-induced cardiomyocyte hypertrophy. The encoded protein has also been found to act as a transcriptional coactivator by interconnecting the general transcription factor TATA element-binding protein (TBP) and gene-specific activators. Alternate splicing results in multiple transcript variants.

*To place an order, please [Click HERE](#).*