

## KAT7

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC IP

**Recommended Dilution:**WB 1:200 - 1:2000 IHC 1:20 - 1:200 IP 1:20 - 1:50

**Calculated MW:**71kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic Peptide of human KAT7

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

f

**Synonym:**

HBO1; HBOA; MYST2; ZC2HC7;

**To place an order, please [Click HERE](#).**

**Catalog #:**A2297

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**11143

**Isotype:**IgG

**Swiss Prot:**O95251

**Purity:**Affinity purification

For research use only.