

KDM3B

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:192kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human KDM3B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

5qNCA; NET22; C5orf7; JMJD1B;

To place an order, please [Click HERE](#).

Catalog #:A2312

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:51780

Isotype:IgG

Swiss Prot:Q7LBC6

Purity:Affinity purification

For research use only.