

SMYD2

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:50kDa

Observed MW:Refer to figures

Immunogen:

A synthetic peptide of human SMYD2

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

KMT3C; HSKM-B; ZMYND14;

Background:

SET domain-containing proteins, such as SMYD2, catalyze lysine methylation (Brown et al., 2006)

To place an order, please [Click HERE](#).

Catalog #:A2327

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:56950

Isotype:IgG

Swiss Prot:Q9NRG4

Purity:Affinity purification

For research use only.