

## SMYD3

---

**Reactivity:**Human

**Tested applications:**WB FC

**Recommended Dilution:**WB 1:200 - 1:2000 FC 1:20 - 1:50

**Calculated MW:**42kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic peptide of human SMYD3

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

f

**Synonym:**

KMT3E; ZMYND1; ZNFN3A1; bA74P14.1;

**Catalog #:**A2329

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**64754

**Isotype:**IgG

**Swiss Prot:**Q9H7B4

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes a histone methyltransferase which functions in RNA polymerase II complexes by an interaction with a specific RNA helicase. Multiple transcript variants encoding different isoforms have been found for this gene.

*To place an order, please [Click HERE](#).*