

MYSM1

Reactivity:Human

Tested applications:WB IHC IP

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:20 - 1:200 IP 1:20 - 1:50

Calculated MW:95kDa

Observed MW:Refer to figures

Immunogen:

A synthetic peptide of human MYSM1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

2ADUB; 2A-DUB;

To place an order, please [Click HERE](#).

Catalog #:A2336

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:114803

Isotype:IgG

Swiss Prot:Q5VVJ2

Purity:Affinity purification

For research use only.