

## TRIB3

---

**Reactivity:**Human Mouse

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:1000

**Calculated MW:**40kDa

**Observed MW:**Refer to Figures

**Immunogen:**

A synthetic peptide of human TRIB3

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

b

**Synonym:**

C20orf97; NIPK; SINK; SKIP3; TRB3;

**Catalog #:**A2346

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**57761

**Isotype:**IgG

**Swiss Prot:**Q96RU7

**Purity:**Affinity purification

For research use only.

**Background:**

Tribbles homolog 3 (TRIB3) is a pseudokinase involved in the regulation of several signaling pathways involved in cell survival and/or cell stress. Whereas mRNA levels of TRIB3 are related with a poor prognosis, TRIB3 protein is associated with a good prognosis in human breast cancer patients, possibly due to the fact that TRIB3 is involved in hypoxia tolerance (PMID:21704407). This antibody is specific to TRIB3.

**To place an order, please [Click HERE](#).**