

Activated NOTCH1

Reactivity:Human Mouse Rat

Tested applications:WB IF FC

Recommended Dilution:WB 1:500 - 1:1000 IF 1:50 - 1:200 FC 1:50 - 1:200

Calculated MW:272kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human Activated NOTCH1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

it

Synonym:

hN1; TAN1;

Catalog #:A2349

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:4851

Isotype:IgG

Swiss Prot:P46531

Purity:Affinity purification

For research use only.

Background:

This gene encodes a member of the Notch family. Members of this Type 1 transmembrane protein family share structural characteristics including an extracellular domain consisting of multiple epidermal growth factor-like (EGF) repeats, and an intracellular domain consisting of multiple, different domain types. Notch family members play a role in a variety of developmental processes by controlling cell fate decisions. The Notch signaling network is an evolutionarily conserved intercellular signaling pathway which regulates interactions between physically adjacent cells. In Drosophila, notch interaction with its cell-bound ligands (delta, serrate) establishes an intercellular signaling pathway that plays a key role in development. Homologues of the notch-ligands have also been identified in human, but precise interactions between these ligands and the human notch homologues remain to be determined. This protein is cleaved in the trans-Golgi network, and presented on the cell surface as a heterodimer. This protein functions as a receptor for membrane bound ligands, and may play multiple roles during development.

To place an order, please [Click HERE](#).