

Pan Acetyl-Lysine rabbit pAb

Reactivity:ALL

Tested applications:WB IP

Recommended Dilution:WB 1:500 - 1:1000 IP 1:20 - 1:50

Observed MW:Refer to Figures

Immunogen:

Acetylated BSA

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A2391

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Isotype:IgG

Purity:Affinity purification

For research use only.

Background:

The rabbit-derived pan anti-acetyllysine polyclonal antibody is purified with acetyl lysine agarose affinity chromatography. It specifically recognizes proteins with acetyllysine residues but not the propionyllysine or butyryllysine residues with structural similarity. Except the common applications such as Western blotting, the antibody is capable of capturing lysine acetylated peptides from protease-digesting proteins both in prokaryotic and eukaryotic cells in systematically proteomic screening.

To place an order, please [Click HERE](#).