

## E2F1

---

**Reactivity:**Human Rat

**Tested applications:**WB IF IP

**Recommended Dilution:**WB 1:200 - 1000 IF 1:50 - 1:100 IP 1:20 - 1:50

**Calculated MW:**47kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human E2F1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

*To place an order, please [Click HERE](#).*

**Catalog #:**A2442

**Antibody Type:**

Monoclonal Antibody

**Species:**Mouse

**Gene ID:**1869

**Isotype:**IgG

**Swiss Prot:**Q01094

**Purity:**Affinity purification

For research use only.