

ATG5

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50- 1:200

Calculated MW:32kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human ATG5

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ASP; APG5; APG5L; hAPG5; APG5-LIKE

To place an order, please [Click HERE](#).

Catalog #:A2468

Antibody Type:

Monoclonal Antibody

Species:Mouse

Gene ID:9474

Isotype:IgG

Swiss Prot:Q9H1Y0

Purity:Affinity purification

For research use only.