

RCN2

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF IP

Recommended Dilution:WB 1:200 - 1:500 IHC 1:50 - 1:100 IF 1:20 - 1:50 IP 1:20 - 1:50

Calculated MW:37kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human RCN2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

E6BP; ERC-55; ERC55; TCBP49;

Catalog #:A2484

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:5955

Isotype:IgG

Swiss Prot:Q14257

Purity:Affinity purification

For research use only.

Background:

Reticulocalbin 2 (RCN2) is a Ca(2+)-binding protein that is specifically localized in the endoplasmic reticulum (ER). It is also termed ERC-55, i.e. ER calcium-binding protein of 55 kDa. It is a single copy gene and is encoded by an approximately 1900-base mRNA, which shows a ubiquitous expression pattern. The protein comprises an amino-terminal signal sequence followed by six copies of the EF-hand Ca²⁺ binding motif.

To place an order, please [Click HERE](#).