

EFNB1

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:38kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human EFNB1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CFND; CFNS; EFL3; EPLG2; Elk-L; LERK2;

Background:

The protein encoded by this gene is a type I membrane protein and a ligand of Eph-related receptor tyrosine kinases. It may play a role in cell adhesion and function in the development or maintenance of the nervous system.

To place an order, please [Click HERE](#).

Catalog #:A2518

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1947

Isotype:IgG

Swiss Prot:P98172

Purity:Affinity purification

For research use only.