

SNCG

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:13kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human SNCG

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SR; BCSG1;

Catalog #:A2524

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:6623

Isotype:IgG

Swiss Prot:O76070

Purity:Affinity purification

For research use only.

Background:

This gene encodes a member of the synuclein family of proteins which are believed to be involved in the pathogenesis of neurodegenerative diseases. Mutations in this gene have also been associated with breast tumor development.

To place an order, please [Click HERE](#).