

HMMR

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:84kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human HMMR

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CD168; IHABP; RHAMM:

Catalog #:A2540

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:3161

Isotype:IgG

Swiss Prot:O75330

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is involved in cell motility. It is expressed in breast tissue and together with other proteins, it forms a complex with BRCA1 and BRCA2, thus is potentially associated with higher risk of breast cancer. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene.

To place an order, please [Click HERE](#).