

## ITGB3

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:**87kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human ITGB3

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

q

**Synonym:**

GT; CD61; GP3A; BDPLT2; GPIIIa; BDPLT16;

**Catalog #:**A2542

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**3690

**Isotype:**IgG

**Swiss Prot:**P05106

**Purity:**Affinity purification

For research use only.

**Background:**

The ITGB3 protein product is the integrin beta chain beta 3. Integrins are integral cell-surface proteins composed of an alpha chain and a beta chain. A given chain may combine with multiple partners resulting in different integrins. Integrin beta 3 is found along with the alpha IIb chain in platelets. Integrins are known to participate in cell adhesion as well as cell-surface mediated signalling.

*To place an order, please [Click HERE](#).*