

GNB2L1

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:35kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human GNB2L1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

jq

Synonym:

H12.3; HLC-7; PIG21; RACK1; Gnb2-rs1;

To place an order, please [Click HERE](#).

Catalog #:A2560

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:10399

Isotype:IgG

Swiss Prot:P63244

Purity:Affinity purification

For research use only.