

SKP1

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW:18kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human SKP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

1 mg/ml

Synonym:

OCP2; p19A; EMC19; SKP1A; OCP-II; TCEB1L;

Catalog #:A2566

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:6500

Isotype:IgG

Swiss Prot:P63208

Purity:Affinity purification

For research use only.

Background:

This gene encodes a component of SCF complexes, which are composed of this protein, cullin 1, a ring-box protein, and one member of the F-box family of proteins. This protein binds directly to the F-box motif found in F-box proteins. SCF complexes are involved in the regulated ubiquitination of specific protein substrates, which targets them for degradation by the proteasome. Specific F-box proteins recognize different target protein(s), and many specific SCF substrates have been identified including regulators of cell cycle progression and development. Studies have also characterized the protein as an RNA polymerase II elongation factor. Alternative splicing of this gene results in two transcript variants. A related pseudogene has been identified on chromosome 7.

To place an order, please [Click HERE](#).