

LAMP1

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:44kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human LAMP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

c

Synonym:

LAMPA; CD107a; LGP120;

Catalog #:A2582

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:3916

Isotype:IgG

Swiss Prot:P11279

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of a family of membrane glycoproteins. This glycoprotein provides selectins with carbohydrate ligands. It may also play a role in tumor cell metastasis.

To place an order, please [Click HERE](#).