

RHOT2

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:68kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human RHOT2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RASL; ARHT2; MIRO2; MIRO-2; C16orf39;

Background:

This gene encodes a member of the Rho family of GTPases. The encoded protein is localized to the outer mitochondrial membrane and plays a role in mitochondrial trafficking and fusion-fission dynamics.

To place an order, please [Click HERE](#).

Catalog #:A2597

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:89941

Isotype:IgG

Swiss Prot:Q8IXI1

Purity:Affinity purification

For research use only.