

OCLN

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 65kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human OCLN

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

OCLN; BLCPMG; Occludin;

Catalog #: A2601

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 100506658

Isotype: IgG

Swiss Prot: Q16625

Purity: Affinity purification

For research use only.

Background:

This gene encodes an integral membrane protein which is located at tight junctions. This protein maybe involved in the formation and maintenance of the tight junction. The possibility of several alternatively spliced products has been suggested but the full nature of these products has not been described. [provided by RefSeq]

To place an order, please [Click HERE](#).