

## LAMTOR5

---

**Reactivity:**Human Mouse

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**10kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant Protein of human LAMTOR5

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

LAMTOR5;HBXIP;XIP;Ragulator complex protein LAMTOR5;

**Background:**

HBXIP, also named as XIP, belongs to the HBXIP family. It is originally identified by its interaction with the C-terminus of hepatitis B virus (HBV) X protein (HBx). HBXIP is a regulator of centrosome dynamics and cytokinesis. HBXIP is able to up-regulate S100A4 through activating STAT4 and inducing DNA methylation of PTEN. HBXIP is a cellular 18kd protein and cytoplasm stain.

*To place an order, please [Click HERE](#).*

**Catalog #:**A2606

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**10542

**Isotype:**IgG

**Swiss Prot:**O43504

**Purity:**Affinity purification

For research use only.