

UPK3B

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:25 - 1:50

Calculated MW:34kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human UPK3B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

P35; UP3B; UPIIIB;

Background:

UPK3B is a minor component of the apical plaques of mammalian urothelium that binds and dimerizes with uroplakin-1b (UPK1B; MIM 602380), one of the major conserved urothelium membrane proteins. The other major conserved integral membrane proteins of urothelial plaques are UPK1A (MIM 611557), UPK2 (MIM 611558), and UPK3A (MIM 611559) (Deng et al., 2002

To place an order, please [Click HERE](#).

Catalog #:A2629

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:80761

Isotype:IgG

Swiss Prot:Q9BT76

Purity:Affinity purification

For research use only.