

AGR2

Reactivity:Human Mouse

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:10 - 1:50 IF 1:20- 1:100

Calculated MW:20kDa

Observed MW:Refer to figures

Immunogen:

A synthetic peptide of human AGR2

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

AG2; GOB-4; HAG-2; XAG-2; PDIA17; HEL-S-116;

To place an order, please [Click HERE](#).

Catalog #:A2632

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:10551

Isotype:IgG

Swiss Prot:O95994

Purity:Affinity purification

For research use only.