

## APIP

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:200 - 1:1000 IHC 1:25 - 1:100

**Calculated MW:**27kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human APIP

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

APIP2; CGI29; hAPIP; CGI-29; MMRP19;

**Background:**

APIP is an APAF1 (MIM 602233)-interacting protein that acts as a negative regulator of ischemic/hypoxic injury (Cho et al., 2004

*To place an order, please [Click HERE](#).*

**Catalog #:**A2644

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**51074

**Isotype:**IgG

**Swiss Prot:**Q96GX9

**Purity:**Affinity purification

For research use only.