

KRT31

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:200 - 1:1000

Calculated MW:47kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human KRT31

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HA1; Ha-1; hHa1; KRTHA1;

Background:

The protein encoded by this gene is a member of the keratin gene family. As a type I hair keratin, it is an acidic protein which heterodimerizes with type II keratins to form hair and nails. The type I hair keratins are clustered in a region of chromosome 17q12-q21 and have the same direction of transcription.

To place an order, please [Click HERE](#).

Catalog #:A2663

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:3881

Isotype:IgG

Swiss Prot:Q15323

Purity:Affinity purification

For research use only.