

## BTG3

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:200 - 1:1000 IHC 1:25 - 1:100

**Calculated MW:**29kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human BTG3

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

ANA; TOB5; TOFA; TOB55;

**Background:**

The protein encoded by this gene is a member of the BTG/Tob family. This family has structurally related proteins that appear to have antiproliferative properties. This encoded protein might play a role in neurogenesis in the central nervous system. Two transcript variants encoding different isoforms have been found for this gene.

*To place an order, please [Click HERE](#).*

**Catalog #:**A2669

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**10950

**Isotype:**IgG

**Swiss Prot:**Q14201

**Purity:**Affinity purification

For research use only.