

CAMK1G

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:53kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human CAMK1G

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

VWS1; CLICK3; CLICKIII; dJ272L16.1;

Background:

This gene encodes a protein similar to calcium/calmodulin dependent protein kinase, however, its exact function is not known.

To place an order, please [Click HERE](#).

Catalog #:A2674

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:57172

Isotype:IgG

Swiss Prot:Q96NX5

Purity:Affinity purification

For research use only.