

## CAP2

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:200 - 1:1000 IHC 1:50 - 1:200

**Calculated MW:**53kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human CAP2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Catalog #:**A2676

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**10486

**Isotype:**IgG

**Swiss Prot:**P40123

**Purity:**Affinity purification

For research use only.

**Background:**

This gene was identified by its similarity to the gene for human adenylyl cyclase-associated protein. The function of the protein encoded by this gene is unknown. However, the protein appears to be able to interact with adenylyl cyclase-associated protein and actin.

**To place an order, please [Click HERE](#).**