

## CDK20

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:1000 IHC 1:50 - 1:100

**Calculated MW:**31kDa

**Observed MW:**Refer to Figures

**Immunogen:**

A synthetic peptide of human CDK20

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

P42; CCRK; CDCH; PNQALRE;

**Background:**

The protein encoded by this gene contains a kinase domain most closely related to the cyclin-dependent protein kinases. The encoded kinase may activate cyclin-dependent kinase 2 and is involved in cell growth. Alternatively spliced transcript variants encoding distinct isoforms have been reported.

*To place an order, please [Click HERE](#).*

**Catalog #:**A2686

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**23552

**Isotype:**IgG

**Swiss Prot:**Q8IZL9

**Purity:**Affinity purification

For research use only.