

CNGA2

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:10 - 1:50

Calculated MW:76kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human CNGA2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CNCA; CNG2; CNCA1; OCNC1; OCNCa; OCNALPHA;

Background:

The protein encoded by this gene represents the alpha subunit of a cyclic nucleotide-gated olfactory channel. The encoded protein contains a carboxy-terminal leucine zipper that mediates channel formation.

To place an order, please [Click HERE](#).

Catalog #:A2698

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1260

Isotype:IgG

Swiss Prot:Q16280

Purity:Affinity purification

For research use only.