

EMC8

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:24kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human EMC8

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NOC4; COX4NB; C16orf2; C16orf4; FAM158B;

To place an order, please [Click HERE](#).

Catalog #:A2701

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:10328

Isotype:IgG

Swiss Prot:O43402

Purity:Affinity purification

For research use only.