

CDX4

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:30kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human CDX4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A2704

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1046

Isotype:IgG

Swiss Prot:O14627

Purity:Affinity purification

For research use only.

Background:

This gene encodes a member of a small subfamily of homeobox containing transcription factors. The encoded protein may regulate homeobox gene expression during anteroposterior patterning and hematopoiesis.

To place an order, please [Click HERE](#).