

IL18R1

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:1000

Calculated MW:62kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human IL18R1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CD218a; IL18RA; IL1RRP; CDw218a; IL-1Rrp;

Background:

The protein encoded by this gene is a cytokine receptor that belongs to the interleukin 1 receptor family. This receptor specifically binds interleukin 18 (IL18), and is essential for IL18 mediated signal transduction. IFN-alpha and IL12 are reported to induce the expression of this receptor in NK and T cells. This gene along with four other members of the interleukin 1 receptor family, including IL1R2, IL1R1, ILRL2 (IL-1Rrp2), and IL1RL1 (T1/ST2), form a gene cluster on chromosome 2q. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).

Catalog #:A2706

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:8809

Isotype:IgG

Swiss Prot:Q13478

Purity:Affinity purification

For research use only.