

CLEC4A

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:1000 IHC 1:50 - 1:200

Calculated MW:28kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human CLEC4A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DCIR; LLIR; DDB27; CLECSF6; HDCGC13P;

Catalog #:A2713

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:50856

Isotype:IgG

Swiss Prot:Q9UMR7

Purity:Affinity purification

For research use only.

Background:

This gene encodes a member of the C-type lectin/C-type lectin-like domain (CTL/CTLD) superfamily. Members of this family share a common protein fold and have diverse functions, such as cell adhesion, cell-cell signalling, glycoprotein turnover, and roles in inflammation and immune response. The encoded type 2 transmembrane protein may play a role in inflammatory and immune response. Multiple transcript variants encoding distinct isoforms have been identified for this gene. This gene is closely linked to other CTL/CTLD superfamily members on chromosome 12p13 in the natural killer gene complex region.

To place an order, please [Click HERE](#).