

DOK3

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:10 - 1:50

Calculated MW:53kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human DOK3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DOKL;

To place an order, please [Click HERE](#).

Catalog #:A2716

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:79930

Isotype:IgG

Swiss Prot:Q7L591

Purity:Affinity purification

For research use only.