

EEF1G

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:100 - 1:200

Calculated MW:50kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human EEF1G

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

EF1G; GIG35;

Catalog #:A2721

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1937

Isotype:IgG

Swiss Prot:P26641

Purity:Affinity purification

For research use only.

Background:

This gene encodes a subunit of the elongation factor-1 complex, which is responsible for the enzymatic delivery of aminoacyl tRNAs to the ribosome. This subunit contains an N-terminal glutathione transferase domain, which may be involved in regulating the assembly of multisubunit complexes containing this elongation factor and aminoacyl-tRNA synthetases.

To place an order, please [Click HERE](#).