

## EGR1

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB

**Recommended Dilution:**WB 1:200 - 1:1000

**Calculated MW:**58kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human EGR1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

TIS8; AT225; GOS30; NGFI-A; ZNF225; KROX-24; ZIF-268;

**Catalog #:**A2722

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**1958

**Isotype:**IgG

**Swiss Prot:**P18146

**Purity:**Affinity purification

For research use only.

**Background:**

The protein encoded by this gene belongs to the EGR family of C2H2-type zinc-finger proteins. It is a nuclear protein and functions as a transcriptional regulator. The products of target genes it activates are required for differentiation and mitogenesis. Studies suggest this is a cancer suppressor gene.

*To place an order, please [Click HERE](#).*